



5067-1522



Declaration of Conformity to Manufacturing Specification

Agilent part number:	5067-1522	Lot number:	DQ19BK10
Product name:	Agilent DNA Chips	Kit expiration date:	19-Mar-2026
Included in:	5067-1504, 5067-1506, 5067-1508		<i>Stewart Knapp</i>
Storage conditions:	Room temperature	Approved by:	Stewart Knapp
Shipping conditions:	Uncontrolled temperature	Date: 31-Mar-2025	Revvity, Inc. Hopkinton, MA 01748

This product was manufactured and verified in accordance with an ISO 9001 certified quality assurance system. The manufacture, inspection, and testing have been performed by trained personnel. The instrumentation and materials used in the qualification process are identified and traceable. All associated test instruments are maintained and calibrated.

During specification verification, the following characteristics were tested on a statistically significant number of chips.

Specification	Method	Result
Channel depth	Profilometer	Pass
Optical inspections	Visual inspections	Pass

For Research Use Only. Not for use in diagnostic procedures.

Agilent Technologies
Deutschland GmbH
Hewlett-Packard-Str. 8
76337 Waldbronn, Germany

Made in USA

For global support: www.agilent.com/en/contact-us/page

Template Number: 5067-1522 DoC
Template Revision: A.00