



5067-1522



## Declaration of Conformity to Manufacturing Specification

|                             |                                 |                             |   |
|-----------------------------|---------------------------------|-----------------------------|---|
| <b>Agilent part number:</b> | 5067-1522                       | <b>Lot number:</b>          | DR17BK10  |
| <b>Product name:</b>        | Agilent DNA Chips               | <b>Kit expiration date:</b> | 17-Apr-2026   |
| <b>Included in:</b>         | 5067-1504, 5067-1506, 5067-1508 |                             | <i>Stewart Knapp</i>                                |
| <b>Storage conditions:</b>  | Room temperature                | <b>Approved by:</b>         | Stewart Knapp                                       |
| <b>Shipping conditions:</b> | Uncontrolled temperature        | <b>Date:</b>                | 07-May-2025<br>Revvity, Inc.<br>Hopkinton, MA 01748 |

This product was manufactured and verified in accordance with an ISO 9001 certified quality assurance system. The manufacture, inspection, and testing have been performed by trained personnel. The instrumentation and materials used in the qualification process are identified and traceable. All associated test instruments are maintained and calibrated.

During specification verification, the following characteristics were tested on a statistically significant number of chips.

| Specification       | Method             | Result |
|---------------------|--------------------|--------|
| Channel depth       | Profilometer       | Pass   |
| Optical inspections | Visual inspections | Pass   |

For Research Use Only. Not for use in diagnostic procedures.

Agilent Technologies  
Deutschland GmbH  
Hewlett-Packard-Str. 8  
76337 Waldbronn, Germany

Made in USA  
For global support: [www.agilent.com/en/contact-us/page](http://www.agilent.com/en/contact-us/page)

Template Number: 5067-1522 DoC  
Template Revision: A.00