

## Declaration of Conformity to Manufacturing Specifications

Product: Agilent Protein Chips

Part No: 5067-1523

Included in Part No: 5067-1515, 5067-1517

Lot No: AT08BP10

---

Note: The chips were prepared using ISO 9001 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

---

### Kit Information

During specification verification, the following characteristics were tested on a statistically significant number of chips:

| Specification       | Method             | Result                                     |
|---------------------|--------------------|--|
| Channel Depth       | Profilometer       | Passed <input checked="" type="checkbox"/> |
| Optical Inspections | Visual Inspections | Passed <input checked="" type="checkbox"/> |
|                     |                    |  |

|  |                |
|--|----------------|
| Kit Expiration Date:                     | 08.JUN.2023    |
| Date:                                    | 28.JUN.2022    |
| Signed:                                  | <b>S.Knapp</b> |
| PerkinElmer, Inc.<br>Hopkinton, MA 01748 |                |

