

## Declaration of Conformity to Manufacturing Specifications

Product: Agilent Protein Chips

Part No: 5067-1523

Included in Part No: 5067-1515, 5067-1517

Lot No: BQ10BP10

---

Note: The chips were prepared using ISO 9001 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

---

### Kit Information

During specification verification, the following characteristics were tested on a statistically significant number of chips:

Specification	Method	Result
Channel Depth	Profilometer	Passed <input checked="" type="checkbox"/>
Optical Inspections	Visual Inspections	Passed <input checked="" type="checkbox"/>

Kit Expiration Date:	10.MAR.2024
Date:	31.MAR.2023
Signed:	<b>K.Garcia</b>
PerkinElmer, Inc. Hopkinton, MA 01748	

