

Declaration of Conformity to Manufacturing Specifications

Product: Agilent Protein Chips

Part No: 5067-1523

Included in Part No: 5067-1515, 5067-1517

Lot No: CP16BP10

Note: The chips were prepared using ISO 9001 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

Kit Information

During specification verification, the following characteristics were tested on a statistically significant number of chips:

Specification	Method	Result
Channel Depth	Profilometer	Passed <input checked="" type="checkbox"/>
Optical Inspections	Visual Inspections	Passed <input checked="" type="checkbox"/>

Kit Expiration Date:	16.FEB.2025
Date:	13.MAR.2024
Signed:	K.Garcia
PerkinElmer, Inc. Hopkinton, MA 01748	

