

Declaration of Conformity to Manufacturing Specifications

Product: Agilent Protein Chips
Included in Part No: 5067-1515, 5067-1517

Part No: 5067-1523
Lot No: ZG23BP10

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. For Research Use Only. Not for use in diagnostic procedures.

Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

Item	Method	Result
Protein Chips		
Channel Depth	Profilometer	Passed <input checked="" type="checkbox"/>
Optical Inspections	Visual Inspections	Passed <input checked="" type="checkbox"/>

Kit Expiration Date	23.JUL.2022
Date: 06.AUG.2021	
Approved by: A.Pfeiffer	
PerkinElmer Hopkinton, MA, USA	

