

Declaration of Conformity to Manufacturing Specifications

Product: Agilent Protein Chips

Part No: 5067-1523

Included in Part No: 5067-1515, 5067-1517

Lot No: ZH25BP10

Note: The chips were prepared using ISO 9001 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

Kit Information

During specification verification, the following characteristics were tested on a statistically significant number of chips:

| Specification | Method | Result |
|---------------------|--------------------|--|
| Channel Depth | Profilometer | Passed <input checked="" type="checkbox"/> |
| Optical Inspections | Visual Inspections | Passed <input checked="" type="checkbox"/> |
| | | |

| | |
|--|-------------|
| Kit Expiration Date: | 25.AUG.2022 |
| Date: 16.SEP.2021 | |
| Signed: S.Knapp | |
| PerkinElmer, Inc. Hopkinton, MA 01748 | |

