

Declaration of Conformity to Manufacturing Specifications

Product: Agilent Protein Chips

Part No: 5067-1523

Included in Part No: 5067-1515, 5067-1517

Lot No: ZL10BP10

Note: The chips were prepared using ISO 9001 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

Kit Information

During specification verification, the following characteristics were tested on a statistically significant number of chips:

| Specification | Method | Result |
|---------------------|--------------------|--|
| Channel Depth | Profilometer | Passed <input checked="" type="checkbox"/> |
| Optical Inspections | Visual Inspections | Passed <input checked="" type="checkbox"/> |
| | | |

| | |
|--|----------------|
| Kit Expiration Date: | 10.NOV.2022 |
| Date: | 06.DEC.2021 |
| Signed: | S.Knapp |
| PerkinElmer, Inc. Hopkinton, MA 01748 | |

