

Declaration of Conformity to Manufacturing Specifications

Product: Agilent RNA Pico Chips
Included in Part No: 5067-1513

Part No: 5067-1524
Lot No: ZG02BK30

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. For Research Use Only. Not for use in diagnostic procedures.

Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

Item	Method	Result
RNA Pico Chips		
Channel Depth	Profilometer	Passed <input checked="" type="checkbox"/>
Optical Inspections	Visual Inspections	Passed <input checked="" type="checkbox"/>

Kit Expiration Date	02.JUL.2022
Date: 21.JUL.2021	
Approved by: A. Pfeiffer	
PerkinElmer Hopkinton, MA, USA	

