

## Declaration of Conformity to Manufacturing Specifications

Product: Agilent High Sensitivity Protein Chips  
Included in Part No: 5067-1575

Part No: 5067-1580  
Lot No: AV12BP20

---

Note: The chips were prepared using ISO 9001 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

---

### Kit Information

During specification verification, the following characteristics were tested on a statistically significant number of chips:

Specification	Method	Result
Channel Depth	Profilometer	Passed <input checked="" type="checkbox"/>
Optical Inspections	Visual Inspections	Passed <input checked="" type="checkbox"/>

Kit Expiration Date:	12.JUL.2023
Date: 20.JUL.2022	
Signed: <b>K. Garcia</b>	
PerkinElmer, Inc. Hopkinton, MA 01748	

