

Declaration of Conformity to Manufacturing Specifications

Product: Agilent High Sensitivity Protein Chips
Included in Part No: 5067-1575

Part No: 5067-1580
Lot No: ZG23BP20

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. For Research Use Only. Not for use in diagnostic procedures.

Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

| Item | Method | Result |
|---------------------------------------|--------------------|--|
| High Sensitivity Protein Chips | | |
| Channel Depth | Profilometer | Passed <input checked="" type="checkbox"/> |
| Optical Inspections | Visual Inspections | Passed <input checked="" type="checkbox"/> |
| | | |

| | |
|------------------------------------|-------------|
| Kit Expiration Date | 23.JUL.2022 |
| Date: 05.AUG.2021 | |
| Approved by: A. Pfeiffer | |
| PerkinElmer Hopkinton, MA, USA | |

