



## CERTIFICATE OF CONFORMANCE

**Agilent Product Name:** A-Line screw cap vials, clear 100/PK

**Agilent Part No:** 5190-9589

**Certificate Date:** January 15<sup>th</sup>, 2017

### Quality Statement

This product is certified and guaranteed to be fully compatible with all corresponding Agilent autosamplers and closures. All lots of this product have been manufactured under a Quality Management System registered to ISO 9001. All glass vials are manufactured in a clean environment, using First Hydrolytic Type 1 borosilicate glass. This product was manufactured to meet a number of critical specifications created by Agilent engineers to ensure a perfect fit with instrumentation.

### Quality Test Results

|                           |  |      |
|---------------------------|--|------|
| Critical Dimensions       | Vial OD                                      | Pass |
|                           | Vial Length                                  | Pass |
|                           | Interior Depth                               | Pass |
|                           | Wall Weight                                  | Pass |
|                           | Thread Major Diameter (GPI- "T" Dim)         | Pass |
|                           | Neck OD/Thread Minor Diameter (GPI- "E" Dim) | Pass |
|                           | Neck ID (GPI- "I" Dim)                       | Pass |
| Optoelectronic Inspection | 100%   | Pass |

Certified by:

Susanne Moyer  
Chemistries and Supplies Division  
Quality Manager