



**Agilent Technologies**

**Certificate of Conformance**

Product: Syringe

Batch No.            LOT25R-052439A

Part No.             5190-1471

Description         Syringe, 2.0 uL, RN 23g, cone tip

Accuracy:             $\pm 2\%$

Reproducibility:     $\pm 1\%$

This syringe has been manufactured to the highest quality using in-process measuring equipment and instrumentation calibrated and traceable to International Standards through the National Research Laboratory of Metrology (NRLM).

Agilent Technologies ensures the displacement is within the stated conformance specifications.