



# Agilent Technologies

## *Certificate of Conformance*

Product: Syringe

***Batch Code***                      ***H09A1535***

***Part Number***                    ***5190-1463***

***Description***                    ***Syringe, 1.0 uL, RN, 26 g, cone tip***

Accuracy:                           $\pm 2\%$

Reproducibility:                  $\pm 1\%$

This syringe has been manufactured to the highest quality using in-process measuring equipment and instrumentation calibrated and traceable to International Standards through the National Research Laboratory of Metrology (NRLM).

Agilent Technologies ensures the displacement is within the stated conformance specifications.