



**Agilent Technologies**

**Certificate of Conformance**

*Product: Syringe*

**Batch Code**    **H04A1325**

**Part**            **5190-1475**

**Description**    **Syringe, 5.0 uL, RN, 23 g, cone tip**

*Accuracy:*             $\pm 2\%$

*Reproducibility:*     $\pm 1\%$

*This syringe has been manufactured to the highest quality using measuring equipment and instrumentation calibrated and traceable to International Standards.*

*Agilent Technologies ensures the displacement is within the stated conformance specifications.*