



# Agilent Technologies

## *Certificate of Conformance*

Product: Syringe

<b>Batch Code</b>	<b><i>H09A1507</i></b>
<b>Part Number</b>	<b><i>5190-1476</i></b>
<b>Description</b>	<b><i>Syringe, 5.0 uL RN 23 g bevel tip</i></b>

Accuracy:	$\pm 2\%$
Reproducibility:	$\pm 1\%$

This syringe has been manufactured to the highest quality using in-process measuring equipment and instrumentation calibrated and traceable to International Standards through the National Research Laboratory of Metrology (NRLM).

Agilent Technologies ensures the displacement is within the stated conformance specifications.