



**Agilent Technologies**

## **Certificate of Conformance**

*Product: Syringe*

**Batch Code**    **H03C9468**

**Part**            **5190-1487**

**Description**    **Syringe, 10.0 uL, FN, 26 g bevel tip, 10/pk**

*Accuracy:*             $\pm 1\%$

*Reproducibility:*     $\pm 1\%$

*This syringe has been manufactured to the highest quality using measuring equipment and instrumentation calibrated and traceable to International Standards.*

*Agilent Technologies ensures the displacement is within the stated conformance specifications.*