



**Agilent Technologies**

## **Certificate of Conformance**

*Product: Syringe*

**Batch Code**    **H05A8826**

**Part**            **5190-1523**

**Description**    **Syringe, 500 uL, PTFE, FN, bevel tip**

*Accuracy:*             $\pm 1\%$

*Reproducibility:*     $\pm 1\%$

*This syringe has been manufactured to the highest quality using measuring equipment and instrumentation calibrated and traceable to International Standards.*

*Agilent Technologies ensures the displacement is within the stated conformance specifications.*