

Agilent Part Number: 5190-9589

Agilent Description: A-Line screw top vial, 2 mL, clear, write-on spot, 100/pk

Certificate Date: February 29, 2024

Quality Statement

This product is verified to be fully compatible with all corresponding Agilent autosamplers and closures. All lots of this product have been manufactured under a Quality Management System registered to ISO 9001. This product was manufactured to meet critical specifications created by Agilent engineers to ensure a compatible fit with instrumentation.

Quality Test Results

Material Type	First Hydrolytic Type 1 Borosilicate Glass	Pass
Critical Dimensions ¹	Vial Outer Diameter	Pass
	Vial Length	Pass
	Interior Depth	Pass
	Wall Weight	Pass
	Thread Major Diameter	Pass
	Neck Outer Diameter / Thread Minor Diameter	Pass
	Neck Inner Diameter	Pass

¹ Confirmed by optoelectronic inspection system

Certified by:



Adam Milecki
Director of Quality Assurance
Chemistries and Supplies Division