

Declaration of Conformity to Manufacturing Specifications

Product: Agilent High Sensitivity DNA Chips
Included in Part No: 5067-4626

Part No: G2938-68000
Lot No: AS25BK50

Note: The chips were prepared using ISO 9001 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

Kit Information

During specification verification, the following characteristics were tested on a statistically significant number of chips:

Specification	Method	Result
Channel Depth	Profilometer	Passed <input checked="" type="checkbox"/>
Optical Inspections	Visual Inspections	Passed <input checked="" type="checkbox"/>

Kit Expiration Date:	25.MAY.2023
Date:	08.JUN.2022
Signed:	S.Knapp
PerkinElmer, Inc. Hopkinton, MA 01748	

