

## Declaration of Conformity to Manufacturing Specifications

Product: Agilent High Sensitivity DNA Chips  
Included in Part No: 5067-4626

Part No: G2938-68000  
Lot No: AX25BK50

---

Note: The chips were prepared using ISO 9001 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

---

### Kit Information

During specification verification, the following characteristics were tested on a statistically significant number of chips:

| Specification       | Method             | Result                                     |
|---------------------|--------------------|--|
| Channel Depth       | Profilometer       | Passed <input checked="" type="checkbox"/> |
| Optical Inspections | Visual Inspections | Passed <input checked="" type="checkbox"/> |
|                     |                    |  |

|  |                 |
|--|-----------------|
| Kit Expiration Date:                     | 25.OCT.2023     |
| Date:                                    | 15.NOV.2022     |
| Signed:                                  | <b>K.Garcia</b> |
| PerkinElmer, Inc.<br>Hopkinton, MA 01748 |                 |

