

Declaration of Conformity to Manufacturing Specifications

Product: Agilent High Sensitivity DNA Chips
Included in Part No: 5067-4626

Part No: G2938-68000
Lot No: BS19BK50

Note: The chips were prepared using ISO 9001 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

Kit Information

During specification verification, the following characteristics were tested on a statistically significant number of chips:

| Specification | Method | Result |
|---------------------|--------------------|--------------------------------------------|
| Channel Depth | Profilometer | Passed <input checked="" type="checkbox"/> |
| Optical Inspections | Visual Inspections | Passed <input checked="" type="checkbox"/> |
| | | |

| | |
|------------------------------------------|----------------|
| Kit Expiration Date: | 19.May.2024 |
| Date: | 05.JUN.2023 |
| Signed: | S.Knapp |
| PerkinElmer, Inc. Hopkinton, MA 01748 | |

