



G2938-68000



## Declaration of Conformity to Manufacturing Specification

<b>Agilent part number:</b>	G2938-68000	<b>Lot number:</b>	CZ23BK50
<b>Product name:</b>	Agilent High Sensitivity DNA Chips	<b>Kit expiration date:</b>	23-Dec-2025
<b>Included in:</b>	5067-4626		<i>Stewart Knapp</i>
<b>Storage conditions:</b>	Room temperature	<b>Approved by:</b>	Stewart Knapp
<b>Shipping conditions:</b>	Uncontrolled temperature	<b>Date:</b> 07-Feb-2025	Revvity, Inc. Hopkinton, MA 01748

This product was manufactured and verified in accordance with an ISO 9001 certified quality assurance system. The manufacture, inspection, and testing have been performed by trained personnel. The instrumentation and materials used in the qualification process are identified and traceable. All associated test instruments are maintained and calibrated.

During specification verification, the following characteristics were tested on a statistically significant number of chips.

Specification	Method	Result
Channel depth	Profilometer	Pass
Optical inspections	Visual inspections	Pass

**For Research Use Only. Not for use in diagnostic procedures.**

Agilent Technologies  
Deutschland GmbH  
Hewlett-Packard-Str. 8  
76337 Waldbronn, Germany

**Made in USA**  
**For global support: [www.agilent.com/en/contact-us/page](http://www.agilent.com/en/contact-us/page)**

Template Number: G2938-68000 DoC  
Template Revision: A.00