

Declaration of Conformity to Manufacturing Specifications

Product: Agilent High Sensitivity DNA Chips
Included in Part No: 5067-4626

Part No: G2938-68000
Lot No: ZJ09BK50

Note: The chips were prepared using ISO 9001 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

Kit Information

During specification verification, the following characteristics were tested on a statistically significant number of chips:

Specification	Method	Result
Channel Depth	Profilometer	Passed <input checked="" type="checkbox"/>
Optical Inspections	Visual Inspections	Passed <input checked="" type="checkbox"/>

Kit Expiration Date:	09.SEP.2022
Date: 21.SEP.2021	
Signed: S.Knapp	
PerkinElmer, Inc. Hopkinton, MA 01748	

