




G2938-80000



## Certificate of Conformity to Manufacturing Specification

<b>Agilent part number:</b> G2938-80000	<b>Lot number:</b> 2403
<b>Product name:</b> Agilent Small RNA Reagents Part I	<b>Kit expiration date:</b> 15-Jan-2025
<b>Included in:</b> 5067-1548 and 5067-1549	
<b>Storage conditions:</b> 2 °C to 8 °C	<b>Approved by:</b> Jurgita Zilinskiene
<b>Shipping conditions:</b> 2 °C to 8 °C	<b>Date:</b> 24-Jan-2024

Thermo Fisher Scientific Baltics  
UAB Vilnius, Lithuania

This product was manufactured and verified in accordance with an ISO 9001 certified quality assurance system. The manufacture, inspection, and testing have been performed by trained personnel. The instrumentation and materials used in the qualification process are identified and traceable. All associated test instruments are maintained and calibrated.

The following characteristics have been tested for conformity to manufacturing specifications.

Product characteristic	Method	Result
<b>Small RNA Gel Matrix</b>		
Conductivity	Electrochemical	Pass
pH	Electrochemical	Pass
No RNase contamination	Electrophoretical	Pass
<b>Small RNA Marker</b>		
Conductivity	Electrochemical	Pass
pH	Electrochemical	Pass
Marker	Electrophoretical	Pass
No RNase contamination	Electrophoretical	Pass
<b>Small RNA Conditioning Solution</b>		
Concentration	Photometrical	Pass
Conductivity	Electrochemical	Pass
pH	Electrochemical	Pass
No RNase contamination	Electrophoretical	Pass
<b>Small RNA Dye Concentrate</b>		
Concentration	Photometrical	Pass

For Research Use Only. Not for use in diagnostic procedures.