




G2938-80000



Certificate of Conformity to Manufacturing Specification

Agilent part number: G2938-80000 **Lot number:** 2528
Product name: Agilent Small RNA Reagents Part I **Kit expiration date:** 07-Jul-2026
Included in: 5067-1548 and 5067-1549
Storage conditions: 2 °C to 8 °C **Approved by:**  Jurgita Zilinskiene
Shipping conditions: 2 °C to 8 °C **Date:** 23-Jul-2025 Thermo Fisher Scientific Baltics
UAB Vilnius, Lithuania

This product was manufactured and verified in accordance with an ISO 9001 certified quality assurance system. The manufacture, inspection, and testing have been performed by trained personnel. The instrumentation and materials used in the qualification process are identified and traceable. All associated test instruments are maintained and calibrated.

The following characteristics have been tested for conformity to manufacturing specifications.

Product characteristic	Method	Result
Small RNA Gel Matrix		
Conductivity	Electrochemical	Pass
pH	Electrochemical	Pass
No RNase contamination	Electrophoretical	Pass
Small RNA Marker		
Conductivity	Electrochemical	Pass
pH	Electrochemical	Pass
Marker	Electrophoretical	Pass
No RNase contamination	Electrophoretical	Pass
Small RNA Conditioning Solution		
Concentration	Photometrical	Pass
Conductivity	Electrochemical	Pass
pH	Electrochemical	Pass
No RNase contamination	Electrophoretical	Pass
Small RNA Dye Concentrate		
Concentration	Photometrical	Pass

For Research Use Only. Not for use in diagnostic procedures.