



G2938-80024



## Certificate of Conformity to Manufacturing Specification

<b>Agilent part number:</b>	G2938-80024	<b>Lot number:</b>	2424
<b>Product name:</b>	Agilent RNA 6000 Pico Reagents Part I	<b>Kit expiration date:</b>	10-Jun-2025
<b>Included in:</b>	5067-1513 and 5067-1514		
<b>Storage conditions:</b>	2 °C to 8 °C	<b>Approved by:</b>	Jurgita Zilinskiene
<b>Shipping conditions:</b>	2 °C to 8 °C	<b>Date:</b>	18-Jun-2024
			Thermo Fisher Scientific Baltics UAB Vilnius, Lithuania

This product was manufactured and verified in accordance with an ISO 9001 certified quality assurance system. The manufacture, inspection, and testing have been performed by trained personnel. The instrumentation and materials used in the qualification process are identified and traceable. All associated test instruments are maintained and calibrated.

The following characteristics have been tested for conformity to manufacturing specifications.

Product characteristic	Method	Result
<b>RNA 6000 Pico Gel Matrix</b>		
Conductivity	Electrochemical	Pass
pH	Electrochemical	Pass
No RNase contamination	Electrophoretical	Pass
<b>RNA 6000 Pico Marker</b>		
Conductivity	Electrochemical	Pass
pH	Electrochemical	Pass
Marker	Electrophoretical	Pass
No RNase contamination	Electrophoretical	Pass
<b>RNA 6000 Pico Conditioning Solution</b>		
Concentration	Photometrical	Pass
Conductivity	Electrochemical	Pass
pH	Electrochemical	Pass
No RNase contamination	Electrophoretical	Pass
<b>RNA 6000 Pico Dye Concentrate</b>		
Concentration	Photometrical	Pass

For Research Use Only. Not for use in diagnostic procedures.

Agilent Technologies  
Deutschland GmbH  
Hewlett-Packard-Str. 8  
76337 Waldbronn, Germany

Made in Lithuania

For global support: [www.agilent.com/en/contact-us/page](http://www.agilent.com/en/contact-us/page)

Template Number: G2938-80024 CoC  
Template Revision: A.00