



G2938-80024



Certificate of Conformity to Manufacturing Specification

Agilent part number: G2938-80024	Lot number: 2431
Product name: Agilent RNA 6000 Pico Reagents Part I	Kit expiration date: 29-Jul-2025
Included in: 5067-1513 and 5067-1514	
Storage conditions: 2 °C to 8 °C	Approved by: Jurgita Zilinskiene
Shipping conditions: 2 °C to 8 °C	Date: 06-Aug-2024

Thermo Fisher Scientific Baltics
UAB Vilnius, Lithuania

This product was manufactured and verified in accordance with an ISO 9001 certified quality assurance system. The manufacture, inspection, and testing have been performed by trained personnel. The instrumentation and materials used in the qualification process are identified and traceable. All associated test instruments are maintained and calibrated.

The following characteristics have been tested for conformity to manufacturing specifications.

Product characteristic	Method	Result
RNA 6000 Pico Gel Matrix		
Conductivity	Electrochemical	Pass
pH	Electrochemical	Pass
No RNase contamination	Electrophoretical	Pass
RNA 6000 Pico Marker		
Conductivity	Electrochemical	Pass
pH	Electrochemical	Pass
Marker	Electrophoretical	Pass
No RNase contamination	Electrophoretical	Pass
RNA 6000 Pico Conditioning Solution		
Concentration	Photometrical	Pass
Conductivity	Electrochemical	Pass
pH	Electrochemical	Pass
No RNase contamination	Electrophoretical	Pass
RNA 6000 Pico Dye Concentrate		
Concentration	Photometrical	Pass

For Research Use Only. Not for use in diagnostic procedures.

Agilent Technologies
Deutschland GmbH
Hewlett-Packard-Str. 8
76337 Waldbronn, Germany

Made in Lithuania

For global support: www.agilent.com/en/contact-us/page

Template Number: G2938-80024 CoC
Template Revision: A.00