

Certificate of Conformance

Product: Agilent Protein 80 Reagents Part II
Included in Part No: 5067-1515 and 5067-1516


Part No: G2938-80040
Lot No: 2142

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

Item	Method	Result
Protein 80 Sample Buffer		
Conductivity	Electrochemical	Passed <input checked="" type="checkbox"/>
pH	Electrochemical	Passed <input checked="" type="checkbox"/>
Marker concentration	Electrophoretical	Passed <input checked="" type="checkbox"/>
Protein 80 Ladder		
Resolution	Electrophoretical	Passed <input checked="" type="checkbox"/>

Kit Expiration Date	18.JUL.2022
Date:	28.10.2021
Approved by:	Jurgita Zilinskiene 
Thermo Fisher Scientific Baltics UAB Vilnius, Lithuania	

