

## Certificate of Conformance

Product: Agilent Protein 80 Reagents Part II  
Included in Part No: 5067-1515 and 5067-1516


Part No: G2938-80040  
Lot No: 2219

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

### Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

Item	Method	Result
<b>Protein 80 Sample Buffer</b>		
Conductivity	Electrochemical	Passed <input checked="" type="checkbox"/>
pH	Electrochemical	Passed <input checked="" type="checkbox"/>
Marker concentration	Electrophoretical	Passed <input checked="" type="checkbox"/>
<b>Protein 80 Ladder</b>		
Resolution	Electrophoretical	Passed <input checked="" type="checkbox"/>

<b>Kit Expiration Date</b>	09.FEB.2023
Date:	23.05.2022
Approved by:	Jurgita Zilinskiene 
Thermo Fisher Scientific Baltics UAB Vilnius, Lithuania	

