

Certificate of Conformance

Product: Agilent Protein 230 Reagents Part II
Included in Part No: 5067-1517 and 5067-1518


Part No: G2938-80041
Lot No: 2208

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

| Item | Method | Result |
|----------------------------------|-------------------|--|
| Protein 230 Sample Buffer | | |
| Conductivity | Electrochemical | Passed <input checked="" type="checkbox"/> |
| pH | Electrochemical | Passed <input checked="" type="checkbox"/> |
| Marker concentration | Electrophoretical | Passed <input checked="" type="checkbox"/> |
| Protein 230 Ladder | | |
| Resolution | Electrophoretical | Passed <input checked="" type="checkbox"/> |
| | | |

| | |
|--|---|
| Kit Expiration Date | 21.NOV.2022 |
| Date: | 02.03.2022 |
| Approved by: | Jurgita Zilinskiene  |
| Thermo Fisher Scientific Baltics UAB Vilnius, Lithuania | |

