

Certificate of Conformance

Product: Agilent Protein 230 Reagents Part II
Included in Part No: 5067-1517 and 5067-1518


Part No: G2938-80041
Lot No: 2235

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

Item	Method	Result
Protein 230 Sample Buffer		
Conductivity	Electrochemical	Passed <input checked="" type="checkbox"/>
pH	Electrochemical	Passed <input checked="" type="checkbox"/>
Marker concentration	Electrophoretical	Passed <input checked="" type="checkbox"/>
Protein 230 Ladder		
Resolution	Electrophoretical	Passed <input checked="" type="checkbox"/>

Kit Expiration Date	29.MAY.2023
Date:	10.08.2022
Approved by:	Jurgita Zilinskiene 
Manufactured for Agilent Technologies	

