




G2938-80041



Certificate of Conformity to Manufacturing Specification

Agilent part number:	G2938-80041	Lot number:	2407
Product name:	Agilent Protein 230 Reagents Part II	Kit expiration date:	12-Nov-2024
Included in:	5067-1517 and 5067-1518		
Storage conditions:	-28 °C to -13 °C	Approved by:	 Jurgita Zilinskiene
Shipping conditions:	2 °C to 8 °C	Date:	27-Feb-2024
			Thermo Fisher Scientific Baltics UAB Vilnius, Lithuania

This product was manufactured and verified in accordance with an ISO 9001 certified quality assurance system. The manufacture, inspection, and testing have been performed by trained personnel. The instrumentation and materials used in the qualification process are identified and traceable. All associated test instruments are maintained and calibrated.

The following characteristics have been tested for conformity to manufacturing specifications.

Product characteristic	Method	Result
Protein 230 Sample Buffer		
Conductivity	Electrochemical	Pass
pH	Electrochemical	Pass
Marker	Electrophoretical	Pass
Protein 230 Ladder		
Resolution	Electrophoretical	Pass

For Research Use Only. Not for use in diagnostic procedures.

Agilent Technologies
Deutschland GmbH
Hewlett-Packard-Str. 8
76337 Waldbronn, Germany

Made in Lithuania

For global support: www.agilent.com/en/contact-us/page

Template Number: G2938-80041 CoC
Template Revision: A.00