




G2938-80041



## Certificate of Conformity to Manufacturing Specification

<b>Agilent part number:</b>	G2938-80041	<b>Lot number:</b>	2446
<b>Product name:</b>	Agilent Protein 230 Reagents Part II	<b>Kit expiration date:</b>	11-Aug-2025
<b>Included in:</b>	5067-1517 and 5067-1518		
<b>Storage conditions:</b>	-28 °C to -13 °C	<b>Approved by:</b>	Jurgita Zilinskiene
<b>Shipping conditions:</b>	2 °C to 8 °C	<b>Date:</b>	20-Nov-2024

Thermo Fisher Scientific Baltics  
UAB Vilnius, Lithuania

This product was manufactured and verified in accordance with an ISO 9001 certified quality assurance system. The manufacture, inspection, and testing have been performed by trained personnel. The instrumentation and materials used in the qualification process are identified and traceable. All associated test instruments are maintained and calibrated.

The following characteristics have been tested for conformity to manufacturing specifications.

Product characteristic	Method	Result
<b>Protein 230 Sample Buffer</b>		
Conductivity	Electrochemical	Pass
pH	Electrochemical	Pass
Marker	Electrophoretical	Pass
<b>Protein 230 Ladder</b>		
Resolution	Electrophoretical	Pass

For Research Use Only. Not for use in diagnostic procedures.

Agilent Technologies  
Deutschland GmbH  
Hewlett-Packard-Str. 8  
76337 Waldbronn, Germany

Made in Lithuania  
For global support: [www.agilent.com/en/contact-us/page](http://www.agilent.com/en/contact-us/page)

Template Number: G2938-80041 CoC  
Template Revision: A.00