

Certificate of Conformance

Product: Agilent High Sensitivity Protein 250 Reagents
Included in Part No: 5067-1575; 5067-1576

Part No: G2938-85000
Lot No: 2207

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

| Item | Method | Result |
|-------------------------|-------------------|--|
| Gel Matrix | | |
| Conductivity | Electrochemical | Passed <input checked="" type="checkbox"/> |
| pH | Electrochemical | Passed <input checked="" type="checkbox"/> |
| Sample Buffer | | |
| Conductivity | Electrochemical | Passed <input checked="" type="checkbox"/> |
| pH | Electrochemical | Passed <input checked="" type="checkbox"/> |
| Marker concentration | Electrophoretical | Passed <input checked="" type="checkbox"/> |
| Destain Solution | | |
| Conductivity | Electrochemical | Passed <input checked="" type="checkbox"/> |
| pH | Electrochemical | Passed <input checked="" type="checkbox"/> |

| | |
|---------------------------------------|---|
| Kit Expiration Date | 14.FEB.2023 |
| Date: | 04.05.2022 |
| Approved by: | Jurgita Zilinskiene  |
| Manufactured for Agilent Technologies | |

