

Certificate of Conformance

Product: Agilent High Sensitivity Protein 250 Reagents
Included in Part No: 5067-1575; 5067-1576


Part No: G2938-85000
Lot No: 2232

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

Item	Method	Result
Gel Matrix		
Conductivity	Electrochemical	Passed <input checked="" type="checkbox"/>
pH	Electrochemical	Passed <input checked="" type="checkbox"/>
Sample Buffer		
Conductivity	Electrochemical	Passed <input checked="" type="checkbox"/>
pH	Electrochemical	Passed <input checked="" type="checkbox"/>
Marker concentration	Electrophoretical	Passed <input checked="" type="checkbox"/>
Destain Solution		
Conductivity	Electrochemical	Passed <input checked="" type="checkbox"/>
pH	Electrochemical	Passed <input checked="" type="checkbox"/>

Kit Expiration Date	08.AUG.2023
Date:	25.08.2022
Approved by:	Jurgita Zilinskiene 
Thermo Fisher Scientific Baltics UAB Vilnius, Lithuania	

