

## Certificate of Conformance

Product: Agilent High Sensitivity Protein 250 Reagents  
Included in Part No: 5067-1575; 5067-1576


Part No: G2938-85000  
Lot No: 2247

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

### Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

Item	Method	Result
<b>Gel Matrix</b>		
Conductivity	Electrochemical	Passed <input checked="" type="checkbox"/>
pH	Electrochemical	Passed <input checked="" type="checkbox"/>
<b>Sample Buffer</b>		
Conductivity	Electrochemical	Passed <input checked="" type="checkbox"/>
pH	Electrochemical	Passed <input checked="" type="checkbox"/>
Marker concentration	Electrophoretical	Passed <input checked="" type="checkbox"/>
<b>Destain Solution</b>		
Conductivity	Electrochemical	Passed <input checked="" type="checkbox"/>
pH	Electrochemical	Passed <input checked="" type="checkbox"/>

<b>Kit Expiration Date</b>	21.NOV.2023
Date:	12.12.2022
Approved by:	Jurgita Zilinskiene 
Thermo Fisher Scientific Baltics UAB Vilnius, Lithuania	

