

Certificate of Conformance

Product: Agilent High Sensitivity Protein 250 Reagents
Included in Part No: 5067-1575; 5067-1576

Part No: G2938-85000
Lot No: 2252

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

Item	Method	Result
Gel Matrix		
Conductivity	Electrochemical	Passed <input checked="" type="checkbox"/>
pH	Electrochemical	Passed <input checked="" type="checkbox"/>
Sample Buffer		
Conductivity	Electrochemical	Passed <input checked="" type="checkbox"/>
pH	Electrochemical	Passed <input checked="" type="checkbox"/>
Marker concentration	Electrophoretical	Passed <input checked="" type="checkbox"/>
Destain Solution		
Conductivity	Electrochemical	Passed <input checked="" type="checkbox"/>
pH	Electrochemical	Passed <input checked="" type="checkbox"/>

Kit Expiration Date	26.DEC.2023
Date:	09.01.2023
Approved by:	Jurgita Zilinskiene 
Manufactured for Agilent Technologies	

