



G2938-85000



Certificate of Conformity to Manufacturing Specification

Agilent part number: G2938-85000**Lot number:** 2615**Product name:** Agilent High Sensitivity Protein 250**Kit expiration date:** 06-Apr-2027

Reagents

Included in: 5067-1575 and 5067-1576**Approved by:**
Jurgita Zilinskiene**Storage conditions:** -28 °C to -13 °C**Date:** 28-Apr-2026Thermo Fisher Scientific Baltics
UAB Vilnius, Lithuania**Shipping conditions:** below -13 °C

This product was manufactured and verified in accordance with an ISO 9001 certified quality assurance system. The manufacture, inspection, and testing have been performed by trained personnel. The instrumentation and materials used in the qualification process are identified and traceable. All associated test instruments are maintained and calibrated.

The following characteristics have been tested for conformity to manufacturing specifications.

Product characteristic	Method	Result
Protein 250 Gel Matrix		
Conductivity	Electrochemical	Pass
pH	Electrochemical	Pass
Protein 250 Sample Buffer		
Conductivity	Electrochemical	Pass
pH	Electrochemical	Pass
Marker	Electrophoretical	Pass
Protein 250 Detain Solution		
Conductivity	Electrochemical	Pass
pH	Electrochemical	Pass

For Research Use Only. Not for use in diagnostic procedures.

Agilent Technologies
Deutschland GmbH
Hewlett-Packard-Str. 8
76337 Waldbronn, Germany

Made in Lithuania

For global support: www.agilent.com/en/contact-us/page

Template Number: G2938-85000 CoC
Template Revision: A.00