

Certificate of Conformance

Product: Agilent High Sensitivity Protein 250 Labeling Reagents
Included in Part No: 5067-1575 and 5067-1577

Part No: G2938-85001
Lot No: 2207

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

| Item | Method | Result |
|------------------------|-----------------|--|
| Labeling Buffer | | |
| pH | Electrochemical | Passed <input checked="" type="checkbox"/> |
| | | |

| | |
|---------------------------------------|---|
| Kit Expiration Date | 14.FEB.2023 |
| Date: | 04.05.2022 |
| Approved by: | Jurgita Zilinskiene  |
| Manufactured for Agilent Technologies | |

