

## Certificate of Conformance

Product: Agilent High Sensitivity Protein 250 Labeling Reagents  
Included in Part No: 5067-1575 and 5067-1577


Part No: G2938-85001  
Lot No: 2228

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

### Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

Item	Method	Result
<b>Labeling Buffer</b>		
pH	Electrochemical	Passed <input checked="" type="checkbox"/>

<b>Kit Expiration Date</b>	11.JUL.2023
Date:	26.07.2022
Approved by:	Jurgita Zilinskiene 
Thermo Fisher Scientific Baltics UAB Vilnius, Lithuania	

