

Certificate of Conformance

Product: Agilent High Sensitivity Protein 250 Labeling Reagents
Included in Part No: 5067-1575 and 5067-1577


Part No: G2938-85001
Lot No: 2232

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

Item	Method	Result
Labeling Buffer		
pH	Electrochemical	Passed <input checked="" type="checkbox"/>

Kit Expiration Date	08.AUG.2023
Date:	25.08.2022
Approved by:	Jurgita Zilinskiene 
Thermo Fisher Scientific Baltics UAB Vilnius, Lithuania	

