




G2938-85001



Certificate of Conformity to Manufacturing Specification

Agilent part number:	G2938-85001	Lot number:	2332
Product name:	Agilent High Sensitivity Protein 250 Labeling Reagents	Kit expiration date:	07-Aug-2024
Included in:	5067-1575 and 5067-1577	Approved by:	 Jurgita Zilinskiene
Storage conditions:	-28 °C to -13 °C	Date:	28-Aug-2023
Shipping conditions:	below -13 °C		Thermo Fisher Scientific Baltics UAB Vilnius, Lithuania

This product was manufactured and verified in accordance with an ISO 9001 certified quality assurance system. The manufacture, inspection, and testing have been performed by trained personnel. The instrumentation and materials used in the qualification process are identified and traceable. All associated test instruments are maintained and calibrated.

The following characteristics have been tested for conformity to manufacturing specifications.

Product characteristic	Method	Result
10 X Protein 250 Labeling Buffer		
pH	Electrochemical	Pass

For Research Use Only. Not for use in diagnostic procedures.

Agilent Technologies
Deutschland GmbH
Hewlett-Packard-Str. 8
76337 Waldbronn, Germany

Made in Lithuania

For global support: www.agilent.com/en/contact-us/page

Template Number: G2938-85001 CoC
Template Revision: A.00