




G2938-85001



Certificate of Conformity to Manufacturing Specification

Agilent part number: G2938-85001 **Lot number:** 2549
Product name: Agilent High Sensitivity Protein 250 **Kit expiration date:** 01-Dec-2026
Labeling Reagents
Included in: 5067-1575 and 5067-1577 **Approved by:**  Jurgita Zilinskiene
Storage conditions: -28 °C to -13 °C **Date:** 05-Dec-2025 Thermo Fisher Scientific Baltics
UAB Vilnius, Lithuania
Shipping conditions: below -13 °C

This product was manufactured and verified in accordance with an ISO 9001 certified quality assurance system. The manufacture, inspection, and testing have been performed by trained personnel. The instrumentation and materials used in the qualification process are identified and traceable. All associated test instruments are maintained and calibrated.

The following characteristics have been tested for conformity to manufacturing specifications.

Product characteristic	Method	Result
10 X Protein 250 Labeling Buffer		
pH	Electrochemical	Pass

For Research Use Only. Not for use in diagnostic procedures.

Agilent Technologies
Deutschland GmbH
Hewlett-Packard-Str. 8
76337 Waldbronn, Germany

Made in Lithuania

For global support: www.agilent.com/en/contact-us/page

Template Number: G2938-85001 CoC
Template Revision: A.00