

Certificate of Conformance

Product: Agilent High Sensitivity Protein 250 Labeling Dye
Included in Part No: 5067-1575, 5067-1577


Part No: G2938-85002
Lot No: E21-05138

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

Item	Method	Result
Labeling Dye		
Purity	HPLC	Passed <input checked="" type="checkbox"/>
Dye Content	UV / Vis	Passed <input checked="" type="checkbox"/>

Kit Expiration Date	17.JUN.2022
Date: 17.JUN.2021	
Approved by: Dr. M. Wenzel	
	
Dyomics GmbH Jena / Germany	

