

## Certificate of Conformance

Product: Agilent High Sensitivity Protein 250 Labeling Dye  
Included in Part No: 5067-1575, 5067-1577

Part No: G2938-85002  
Lot No: E21-06109

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

### Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

Item	Method	Result
<b>Labeling Dye</b>		
Purity	HPLC	Passed <input checked="" type="checkbox"/>
Dye Content	UV / Vis	Passed <input checked="" type="checkbox"/>

<b>Kit Expiration Date</b>	07.DEZ.2023
Date: 07.DEZ.2022	
Approved by: Dr. M. Wenzel	
	
Dyomics GmbH Jena / Germany	

