

## Certificate of Conformance

Product: Agilent High Sensitivity Protein 250 Labeling Dye  
Included in Part No: 5067-1575, 5067-1577


Part No: G2938-85002  
Lot No: E21-06191

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

### Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

Item	Method	Result
<b>Labeling Dye</b>		
Purity	HPLC	Passed <input checked="" type="checkbox"/>
Dye Content	UV / Vis	Passed <input checked="" type="checkbox"/>

<b>Kit Expiration Date</b>	12.OCT.2024
Date: 12.OCT.2023	
Approved by: Dr. M. Wenzel	
	
Dyomics GmbH Jena / Germany	

