




G2938-85002



Certificate of Conformity to Manufacturing Specification

Agilent part number:	G2938-85002	Lot number:	E21-07039
Product name:	Agilent High Sensitivity Protein 250 Labeling Dye	Kit expiration date:	06.FEB.2025
Included in:	5067-1575 and 5067-1577	Approved by:	 Dr. M. Wenzel
Storage conditions:	-28 °C to -13 °C	Date:	06.FEB.2024 Dyomics GmbH, Jena Germany
Shipping conditions:	below -13 °C		

This product was manufactured and verified in accordance with an ISO 9001 certified quality assurance system. The manufacture, inspection, and testing have been performed by trained personnel. The instrumentation and materials used in the qualification process are identified and traceable. All associated test instruments are maintained and calibrated.

The following characteristics have been tested for conformity to manufacturing specifications.

Product characteristic	Method	Result
High Sensitivity Protein 250 Labeling Dye		
Purity	High Performance Liquid Chromatography	Pass
Dye Content	UV/VIS Spectroscopy	Pass

For Research Use Only. Not for use in diagnostic procedures.

Agilent Technologies
Deutschland GmbH
Hewlett-Packard-Str. 8
76337 Waldbronn, Germany

Made in Germany

For global support: www.agilent.com/en/contact-us/page

Template Number: G2938-85002 CoC
Template Revision: A.00