

Certificate of Conformance

Product: Agilent High Sensitivity Protein 250 Ladder
Included in Part No: 5067-1575, 5067-1578


Part No: G2938-85003
Lot No: 2139

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

Item	Method	Result
High Sensitivity Protein 250 Ladder		
Concentration	Photometrical	Passed <input checked="" type="checkbox"/>
Functional	Electrophoretical	Passed <input checked="" type="checkbox"/>

Kit Expiration Date	27.SEP.2022
Date:	29.09.2021
Approved by:	Jurgita Zilinskiene 
Thermo Fisher Scientific Baltics UAB Vilnius, Lithuania	

