

Certificate of Conformance

Product: Agilent High Sensitivity Protein 250 Ladder
Included in Part No: 5067-1575, 5067-1578


Part No: G2938-85003
Lot No: 2235

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

| Item | Method | Result |
|--|-------------------|--|
| High Sensitivity Protein 250 Ladder | | |
| Concentration | Photometrical | Passed <input checked="" type="checkbox"/> |
| Functional | Electrophoretical | Passed <input checked="" type="checkbox"/> |
| | | |

| | |
|---------------------------------------|---|
| Kit Expiration Date | 29.AUG.2023 |
| Date: | 14.09.2022 |
| Approved by: | Jurgita Zilinskiene  |
| Manufactured for Agilent Technologies | |

